

Bayer® Glucose Meter Characteristics

	Bayer® CONTOUR®	Bayer® BREEZE® 2
		
Features	<ul style="list-style-type: none"> -Auto Calibration -5 second test time -Underfill detection -Before and after meal markers -2 Settings-Basic(L1) & Advanced (L2) modes Automatic hematocrit correction -Hematocrit range 0-70% -Humidity range 10-93% -Unaffected by maltose and galactose 	<ul style="list-style-type: none"> -Auto Calibration -5 seconds test time -Underfill detection -Operating Temperature 50°-113°F -Hematocrit range 20-55% -Humidity range 10-80% -Discreet pull-push mechanism -No manual cleaning required -Unaffected by maltose and galactose
Blood Sample size	Only 0.6µL	Small sample size of 1µL
Alternate site testing	Approved	Approved
Memory	480 tests with date & times 7-day Hi/Lo test summary Individual HI and Lo ranges 7-,14-, 30- day averages	420 tests with date & times 1,7,14, and 30 day test averages
PC download capability	GLUCOFACTS® Deluxe Diabetes Management Software can be downloaded at http://www.simplewins.com/site/Adults/Monitor/Pages/Bayer-Products/GLUCOFACTS-Deluxe-Software . Download your test results from your CONTOUR® or BREEZE®2. User-friendly charts and graphs to help you find patterns in your blood sugar tests. CONTOUR® USB Meter, the first blood glucose meter with plug and play diabetes management software, available only with a prior authorization.	
Test strips	CONTOUR® individual test strips. Sip-in Sampling® (with capillary action)	BREEZE2 Test Strips (10-test disc)
Use in patients with EXTRANEAL (icodextrin) peritoneal dialysis solution.	BOTH METERS HAVE APPROVAL FOR USE Test results not affected by icodextrin.	
Design	Small meter Large Display and buttons Weight of 1.7 ounces	 Single function buttons, especially for people with limited dexterity. Comfortable and convenient shape
Size (LXWXH) inches	3.0" x 2.2" x 0.75"	4.3" x 2.3" x 1.1"
Batteries	1000 tests/All Models Batteries Free For Life. Call 1- 800-348-8100	
Warranty	5 year warranty	
Customer Support	Call 1-800-348-8100 for 24/7 customer support.	

To obtain a FREE Bayer meter, members or providers can call toll free to 800-348-8100 or visit

<https://offers.bayerdiabetes.com/CouponReg.aspx?product=contour>

Demonstrations of Bayer® products and current product updates are available on-line at www.bayerdiabetes.com.